

# **PLANT GROWTH CHAMBER**

1


RK-340 CH RK-500 CH RK-1000 CH

Perfect solution for botanical and agricultural applications, plant pathology, plant tissue culture studies, seed germination studies, plant genetic research...

Three different sizes (small, medium, large)

- Easy to use fully programmable controller
  - State of the art repeatability Day light simulation control Relative humidity control Temperature control
    - CO2 control (option)







#### **Controller functions:**

#### 24 hour day & night cycle control:

24 hours (one day) can be divided into five phases of day and five phases of night!


Example:

- Phase 1 (day): 07:00 to 10:00, T = 15°C, Rh = 70%, light intensity 1
- Phase 2 (day): 10:00 to 18:00, T = 30°C, Rh = 80%, light intensity 3
- Phase 3 (day): 18:00 to 21:00, T = 23°C, Rh = 75%, light intensity 2


- .....

- Phase 1 (night): 21:00 to 01:00, T = 20°C, Rh = 70%, no lights
- Phase 2 (night): 01:00 to 05:00, T = 18°C, Rh = 73%, no lights
- Phase 3 (night): 05:00 to 07:00, T = 19°C, Rh = 75%, no lights





Auto recovery after power cut-out:





### Technical data:

	RK-340 CH	RK-500 CH	RK-1000 CH
External dimensions (WxHxD) in mm	834 x 1895 x 1135	1600 x 1980 x 700	2040 x 1983 x 886
Internal dimensions (WxDxH) in mm	600 x 830 x 685	680 x 1400 x 520	1045 x 1400 x 700
Volume (L)	~340	~500	~1000
Number of shelves with lights (lights on shelves option)	2	3	3
Number of wire shelves	2	3	3///
Temperature range (°C) (lights on or off)	+5+60	+5+60	+5+60
Rh range (%) (DEW point limits)	1095	1095	1095
Temperature & Rh display resolution	0.1°C / 1%	0.1°C / 1%	0.1°C/1%
Temperature set resolution (°C)	0.1	0.1	0.1
Rh set resolution (%)	1	1//////////////////////////////////////	
Temperature stability (°C) (lights off)	Better than ±0,3°C	Better than ±0,3°C	Better than ±0,3°C
Temperature uniformity (°C) (lights off)	< ± 1 @ +37°C	< ± 1 @ +37°C	< ± 1 @ +37°C
Rh Stability (%) (lights off)	Better than ±3	Better than ±3	Better than ±3
Heating rate	~ 2°C/min	~/2°C/min	~ 2°C/min
Cooling rate	~ 1.5°C/min	~ 1.5°C/min	~ 1.5°C/min
Day light simulation	FLUORA lights (LED option)	FLUORA lights (LED option)	FLUORA lights (LED option)
Illumination intensity	1/3, 2/3, 3/3	1/3, 2/3, 3/3	1/3, 2/3, 3/3
Temperature Control	PID	PUP THPID	PID
Cooling	Air cooled compressor	Air cooled compressor	Air cooled compressor
	aggregate	aggregate	aggregate
Power supply	3x400 V 50/60 Hz (± 10%)	3x400 V 50/60 Hz (± 10%)	3x400 V 50/60 Hz (± 10%)
Wattage (W)	4500	4500	4500
Interface	RS 232 (USB or Ethernet	RS/232 (USB or Ethernet	RS 232 (USB or Ethernet
Maight (kg)	optional)	optional)	optional)
Weight (kg)	318///	469	642



## Ordering information and accessories:

Description	Part no.	
Plant Growth Chamber RK-340 CH	922	
Plant Growth Chamber RK-500 CH	903	
Plant Growth Chamber RK-1000 CH	2638	
Access Port w. Plug Ø 40 mm	2629	
Access Port w. Plug Ø 50 mm	608////	
Access Port w. Plug Ø 90 mm	1731	
Evaluation report	1719///	
Monitoring and data collection software for PC	104	
USB interface	1466	
Ethernet interface	1716	
Automatic water feeding	1744	
Password protection	1718	
CO2 control system	2108	
Shelve wire RK-340 CH	2109	
Shelve wire RK-500 CH	/2110	
Shelve wire RK-1000 CH	/2111	
Shelve with lights RK-340 CH	2112	
Shelve with lights RK-500 CH	2113	
Shelve with lights RK-1000 CH	2114	



