

“WALK IN” PLANT GROWTH CHAMBER

Perfect solution for botanical and agricultural applications, plant pathology, plant tissue culture studies, seed germination studies, plant genetic research...

- Custom sizes according to user requirements
- Easy to use fully programmable controller
 - State of the art repeatability
 - Day light simulation control
 - Relative humidity control
 - Temperature control
 - CO2 control (option)
 - PC compatibility






Controller functions:

24 hour day & night cycle control:


24 hours (one day) can be divided into five phases of day and five phases of night!

Example:


- Phase 1 (day): 07:00 to 10:00, $T = 15^{\circ}\text{C}$, $R_h = 70\%$, light intensity 1
- Phase 2 (day): 10:00 to 18:00, $T = 30^{\circ}\text{C}$, $R_h = 80\%$, light intensity 3
- Phase 3 (day): 18:00 to 21:00, $T = 23^{\circ}\text{C}$, $R_h = 75\%$, light intensity 2
-
- Phase 1 (night): 21:00 to 01:00, $T = 20^{\circ}\text{C}$, $R_h = 70\%$, no lights
- Phase 2 (night): 01:00 to 05:00, $T = 18^{\circ}\text{C}$, $R_h = 73\%$, no lights
- Phase 3 (night): 05:00 to 07:00, $T = 19^{\circ}\text{C}$, $R_h = 75\%$, no lights
-



Auto recovery after power cut-out:



Technical data:

	Walk in Plant Growth Chamber
	
External dimensions (WxHxD) in mm	According to user preferences
Temperature range (°C) (lights on or off)	-20...+60
Rh range (%)	40...95
Temperature & Rh display resolution	0.1°C / 1%
Temperature set resolution (°C)	0.1
Rh set resolution (%)	1
Temperature stability (°C) (lights off)	Better than ±0,3
Temperature uniformity (°C) (lights off)	< ± 1 @ +37°C
Rh Stability (%) (lights off)	Better than ±3
Day light simulation	FLUORA lights (LED as an option)
Illumination intensity	1/3, 2/3, 3/3
Temperature Control	PID
Cooling	Air cooled compressor aggregate
Relative humidity	Steam generator / Ultrasonic humidifier / DEW point
Power supply	3x400 V 50/60 Hz (± 10%)
Access ports (mm)	Ø 40 standard, (Ø 50, Ø 90 as an accessory)
Interface	RS 232 (USB or Ethernet as an optional)

Accessories:

Description	Part no.
Evaluation report	1719
Monitoring and data collection software for PC	104
USB interface	1466
Ethernet interface	1716
Password protection	1718
CO2 control system	2108

Design solutions:

Open room installation:



Incorporation in to a Clean Room:



Outside installation:



Portable:
(Incorporated in to 20ft or 40ft
ship container)

